



To Whom It May Concern:

Please find below, information related to your client's hearing needs.

Venues that your client attends are equipped with a hearing loop system. This system provides the capability for telecoil equipped hearing aids to directly connect with the sound transmission within the facility. For hearing aids to effectively connect with the loop system, certain requirements need to be met. Please find these listed below:

1. The telecoil should be oriented in the vertical direction as worn in or on the ear or in the streamer/gateway device.
2. The hearing aid wearer must be able to turn on and control the telecoil manually.
3. Optimally, a T and an M/T program should be provided. The M output for the M/T program should be 3 – 6 dB lower than the T output.
4. The telecoil performance should be measured via an HFA-SPLIV test such that a 31.6mA/m input to the telecoil provides the same output as the microphone provides with a 60dB SPL input.
5. ***Complete testing and client training with the hearing loop system installed in your clinic.***

A more detailed description of the setup procedure can be found at http://www.audiologyonline.com/articles/article_detail.asp?article_id=2388 - especially points 11 - 16.

Feel free to contact me for clinic hearing loop systems or any additional information that you may require.

Sincerely,

W. Droogendyk

